INSTRUCTIONS

CRYSTAL CLEAR RINGMASTER™

IR OVEN INSTRUCTIONS & NOTE FOR CONVECTION OVEN

TIMES AND TEMPS ARE APPROXIMATE

1. Pre-Heat oven to 365° F
2. Prepare positive model by setting up fixture on distal end (such as TF-65) and use plaster to blend the edges.
3. Perform regular prep on positive model.
4. Place aluminum support ring (RM-AS) in oven supporting it by at least two sides (RM-ST support stand shown.) Make sure the inside diameter of the support ring is open.
5. Place Ringmaster™ sheet on the aluminum support, then place second aluminum support ring (RM-AS) on top, centered with plastic.
6. Heat inside oven for approximately 15 minutes.
7. When plastic droops (1/2-2/3 of model,) open oven and pull the rack out. Grab the Ringmaster™ and frame, take to model, flip like normal blister forming, pull down with a slight twisting motion and apply vacuum.
8. Use a flat head screw driver and place between support ring and plastic, then twist, creating a void. This will allow for easy separation of the support ring and the plastic when the plastic has cooled. Remove the ring.
9. Make trim lines and cut the socket.

Note: When using in a convection oven, flip and/or rotate the plastic as necessary to keep even wall thickness.

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